

Scintillator-Layered Imaging Microscope for Environmental Research

M.F. Kidd, S.R. Elliott, (Los Alamos National Laboratory)


T.C. Onstott, S. Myeni, B. Stackhouse, (Princeton University)

S.M. Pfiffner, T. Vishnivetskaya, A. Layton, (University of Tennessee)

L.G. Whyte, N. Mykytczuk, J. Allan, R.C. Wilhem, (McGill University)

R. Hettich, K. Chourey, T.J. Phelps, (Oak Ridge National Laboratory)

P. Hatcher, (Old Dominion University)



Scintillator-Layered Imaging Microscope for Environmental Research

SLIMER

M.F. Kidd, S.R. Elliott, (Los Alamos National Laboratory)

T.C. Onstott, S. Myeni, B. Stackhouse, (Princeton University)

S.M. Pfiffner, T. Vishnivetskaya, A. Layton, (University of Tennessee)

L.G. Whyte, N. Mykytczuk, J. Allan, R.C. Wilhem, (McGill University)

R. Hettich, K. Chourey, T.J. Phelps, (Oak Ridge National Laboratory)

P. Hatcher, (Old Dominion University)

LA-UR-11-10260

Studying the Permafrost

- * Permafrost underlies 22-24% of Earth's surface.
- * Global climate models project the strongest future warming in Northern Hemisphere high latitudes.
- * Thawing permafrost results in microbial decomposition of previously frozen organic C.

Studying the Permafrost

- * We need to understand the rate at which individuals in the microbial community ingest the organic C.
- * The microbial community will be fed enriched ^{14}C , which decays via beta decay.
- * The microbial community will be divided by type onto a RNA/DNA microarray.
- * Each type will be fluorescently tagged.



SLIMER's Objective

- * To image fluorescent light from tagged biological samples.
- * To detect beta decays from ^{14}C contained within the biological samples.
- * **This has never been done in a single step process.**

Previous System - Isotope Array

- * Used a beta-imager for radioactivity analysis
- * two step process
- * limited resolution

APPLIED AND ENVIRONMENTAL MICROBIOLOGY, Nov. 2003, p. 6875–6887
0099-2240/03/\$08.00+0 DOI: 10.1128/AEM.69.11.6875–6887.2003
Copyright © 2003, American Society for Microbiology. All Rights Reserved.

The Isotope Array, a New Tool That Employs Substrate-Mediated Labeling of rRNA for Determination of Microbial Community Structure and Function

Justyna Adamczyk,¹ Martin Hesselsoe,² Niels Iversen,² Matthias Horn,³ Angelika Lehner,¹ Per Halkjaer Nielsen,² Michael Schloter,⁴ Peter Roslev,^{2*} and Michael Wagner³

Lehrstuhl für Mikrobiologie, Technische Universität München, D-85350 Freising,¹ and GSF-Forschungszentrum für Umwelt und Gesundheit GmbH, D-85764 Neuherberg,⁴ Germany; Section of Environmental Engineering, Aalborg University, DK-9000 Aalborg, Denmark;² and Abteilung für Mikrobielle Ökologie, Institut für Ökologie und Naturschutz, Universität Wien, A-1090 Vienna, Austria³

Received 7 May 2003/Accepted 14 August 2003





FIG. 2. Isotope array experiment with *N. europaea* Nm57. rRNA was extracted from *N. europaea* cultures after incubation with [¹⁴C]bicarbonate. After labeling of the rRNA with Cy5, a single slide carrying three prototype arrays with different spot sizes was hybridized, washed, and scanned for fluorescence and radioactivity. (A) Spot diameter, 125 μm. (B) Spot diameter, 500 μm. (C) Approximate spot diameter, 1 mm. Bars, 500 μm.

SLIMER's Objective



SLIMER's Objective

- * Idea: overlay biological samples with thin scintillator.
 - * electrons from ^{14}C produce light in scintillator
 - * fluorescent light passes through
- * Therefore, need scintillator with:
 - * correct optical properties
 - * ability to maintain location information

CsI(Tl)

- * Thallium-doped cesium iodide (CsI(Tl)) is a commonly used scintillator in nuclear physics
- * Light output is about 65 photons/keV peaked at 565 nm.
- * When grown on a substrate with a structure (i.e. fiber optic plate) the CsI(Tl) forms columns which channel the light
- * Transmits a large range of wavelengths.




Figure 1. SEM of a 140 μm thick micro-columnar CsI(Tl) film fabricated at RMD


CsI(Tl)




45 MICRON

135 MICRON

150 MICRON




Edge of the CsI(Tl)







System Components

- * Purchased Nikon Eclipse Ti-U.
- * Considered transmission of each component:
 - * Objectives >85%
 - * Filter cube >95%
- * Photometrics Evolve- EMCCD
 - * on-board EM gain
 - * cooled as low as -85 °C




The Setup





The Signal

EMCCD




CsI(Tl)

Sample



The Signal

EMCCD




CsI(Tl)




Sample

The Signal

EMCCD




The Signal




Analysis

- * Two potential methods:
 - * frame-by-frame
 - * long-exposure



Analysis - Frame- by-Frame

- * ^{241}Am source
- * 5 MeV alphas make obvious signal



Analysis - Frame- by-Frame

- * Subtract background
- * Gaussian smear





Analysis - Frame- by-Frame

- * Apply threshold




Analysis




Cs-137

- * Electrons emitted in beta decay have endpoint of 513.97 keV (94.70%) with mean energy of 174.32 keV
- * Deposit between 80 to 108 keV in 150 μm of CsI
- * (Recall ^{14}C has endpoint of 156.48 keV and mean of 49.47 keV)




Spectra

- * ^{137}Cs with 250 μm mask.




Spectra

- * Background with 250 μm mask.





Spectra





Long Exposures

- * Can also take long exposures to image shapes.




Long Exposures

- * Can also take long exposures to image shapes.





Long Exposures

- * 250 μm mask with 3 minute exposure.




Sr-90


- * emits a 195.8 keV (mean), 546.0 keV (endpoint) beta to ^{90}Y .
- * 195.8 keV electron should deposit about 109.4 keV in 150 μm of CsI.



59 keV gamma from ^{241}Am




C-14?




45 nCi source

C-14?




C-14?




The background problem

- * As seen, it's very difficult to pull out the ^{14}C signal.
- * We have asked around the lab if anyone has a stronger ^{14}C source for testing.
- * Meanwhile, we continue to investigate the noise in the camera




The background problem

- * Naively, one would expect background to be lower at lower temperatures
- * This is generally true, however EMCCDs have a different feature - clock-induced charge (CIC)
- * At low temperatures and high gain, this is a noticeable effect.




BERT

- * One way to mitigate CIC is to utilize a tool in the microscope software - Background Event Reduction Technology (BERT)
- * BERT sets up a median filter to remove spurious events in real time
 - * each pixel is compared to surrounding pixels
 - * if the pixel is greater than the median value of itself and surrounding pixels by a user-defined threshold, it is replaced by the median value




BERT side-by-side




To Do

- * Detect ^{14}C convincingly
- * Continue to test resolution
- * Continue to improve frame-by-frame analysis
- * Study sample slide - should arrive this week!


Thank you!



^{241}Am
(4.5 nCi)



Background



^{137}Cs
(5 μCi)

Supported by DOE-BER